


Fusion PP+



Suggested fitting range:



Fusion PP+ Body Aid

Technical Description:

The Fusion PP+ is a very small, highly powerful, compact body aid which easily fits in the breast pocket of a shirt.

Fitting Recommendation:

The Fusion PP+ is recommended for mild to profound hearing losses.

Standard Features:

- Small, pager-like housing
- Black cord and receiver
- Battery type AA
- Large volume control dial numbered 1-8
- Reversible clip
- Telecoil
- Direct audio input socket
- Tone control trimmer (NH)
- Automatic gain control (AGC-O)

Available Options:

- Harness for pediatric fitting
- Y-cord for binaural fitting
- Various cord lengths
- Wide band receiver

Warning:

This instrument has an output sound pressure level that may exceed 132 dB. Therefore, special care should be exercised in the selection and fitting of this instrument to avoid impairing the remaining hearing of the patient.

ANSI S3.22-1996

MAXIMUM POWER OUTPUT	Peak	dB OSPL90	140
	HF-average	dB OSPL90	136
GAIN (Input at 50 dB)	Peak	dB	80
	HF-average	dB	76
FREQUENCY RESPONSE		Hz	130-3800
COMPRESSION	Attack time	ms	140
	Release time	ms	600
TOTAL HARMONIC DISTORTION	Maximum	%	<8
EQUIVALENT NOISE LEVEL		dB SPL	<27
INDUCTION COIL SENSITIVITY	10 mA/m	dB SPL	124
BATTERY DRAIN	at RTG., Input 65 dB SPL	mA	7
BATTERY LIFE	Type AA	hrs	280

5010 Cheshire Ln. N., Suite 2
Plymouth, MN 55446
800-876-1141
Fax: 763-553-9129
www.rexton-online.com

REXTON

Fusion PP+ Body Aid

