

INTUIS™



Easy on the ears.

Siemens Hearing Instruments
Better hearing is better living.

SIEMENS

INTUIS Easy to wear.

The new INTUIS™ family of hearing instruments from Siemens sets new standards in hearing ease and comfort by offering you key advancements in performance and reliability at a tremendous value. INTUIS comes in a wide range of models and colors to fit most hearing losses and personal styles.

The attractive directional Behind-the-Ear instruments include power and superpower models. The custom line-up offers almost invisible in-the-ear Micro-CIC™ models with more power than ever before. The INTUIS family even features a very cosmetically appealing open fit BTE. Effective, reliable and simply intuitive – you'll be amazed at just how easy hearing is with INTUIS.



Intuitive – inside and out.

In addition to personal design and comfort, INTUIS integrates select leading technologies in reliability and quality. With proprietary technologies to protect from cerumen (ear wax), sweat and moisture, as well as manufacturing innovations to ensure long-term performance, no other comparable solution offers as much as INTUIS.



INTUIS S Dir



INTUIS SP Dir



INTUIS Dir



INTUIS Life



INTUIS instruments are available in a variety of colors to match virtually every hair color and skin complexion.



INTUIS
ITC/Half-Shell



INTUIS
Mini Canal/
Power MC



INTUIS
Micro-CIC/
Power CIC



INTUIS
ITE/Power ITE

INTUIS Technology at a glance.

Directional microphone system

Enhances speech understanding in noisy places. Available on select models.

Feedback cancellation

Eliminates or reduces feedback (whistling) that can occur with other instruments.

Noise reduction

Effectively reduces steady-state background noises for better comfort.

AutoPhone

Allows you to conveniently and naturally use the phone. Available on select custom models.

Nanocoating

Unique protection for BTE models from moisture and debris.

C-Guard™

Revolutionary wax-protection membrane option for custom instruments.

Hearing instruments may help people hear better, but no hearing instrument can solve every problem or restore hearing back to normal.

Information in this brochure is subject to change without notice.
Copyright ©2007 Siemens Hearing Instruments, Inc. All rights reserved.

Siemens Hearing Instruments, Inc.
P.O. Box 1397
Piscataway, NJ 08855-1397
www.usa.siemens.com/hearing

Siemens Hearing Instruments
A Division of Siemens Canada Limited
320 Pinebush Road,
Cambridge, Ontario N1T 1Z6
www.siemens.ca/hearing

INTUIS is available at: