ACCESSORIES PRICE LIST

ILLUSTRATIONS ARE NOT ACTUAL SIZE

Foam Telephone Pad..\$3.00 Bulb Air Blower\$6.00

Comply Ear Tip'sea. \$4.00 Wax Picks\$2.50

Pouch\$4.00

Soft

Wraps......\$9.95/pkg. 10

BATTERY CHECKER/ HOLDER KEYCHAIN



\$1.50



FOR BEHIND-THE-EAR AIDS



Please specify your Hearing Aid model

\$2.00#13C **Curved Tubing** for Custom

Earmolds

\$2.00 To Remind the Ear \$2.00/Ft. Core & #13T Tubing Straight for Inserts Tubina

MUSHROOM **INSERTS**



TLM-SS ...X-Small TLM-SSmall TLM-M....Medium TLM-LLarge TLM-LL....X-Large

DRI-AID KIT

Moisture is the enemy of hearing aids. Protect against corrosion and deterioration by removing the damaging moisture with Dri-Aid.

\$8.00

Inexpensive and very effective! State choice of pouch or canister.

FOR BODY AIDS



Receivers: Specify Model-\$30.00 **Body Aid**

\$2.00

#TLC

\$2.50Mushroom Insert


Connector (See above)



CORDS

Straight Cord - \$20.00 Y-Cord - \$26.00

Please specify 2 or 3 prong and specify length and your aid model (18", 24", 30", or 36")



TROUBLE WITH WAX GETTING INTO YOUR AID? HERE'S THE ANSWER

Earwax getting inside the aid is by far the #1 cause for repairs of I-T-E hearing aids. Avoid costly repairs and downtime with your instruments. This proven remedy works similar to a "Band-aid" as it covers the receiver outlet hole with a special pass-through filter allowing the sound to pass clearly while keeping out damaging earway or wetness. Each