



Introducing the new Targa *PLUS* and Targa *PRO* from Rexton

Today's hearing technology has brought forth a new era of improved hearing. When compared to traditional hearing instruments, the new Targa series offers more precise instrument fittings and produces excellent sound quality. These advanced hearing instruments are an excellent choice for individuals looking for a better quality of life.

The completely redesigned Targa family of hearing instruments offers unique hearing solutions that incorporate more advanced digital features and user-driven options than ever before. All Targa instruments are economically priced to ensure that budget-conscious consumers can receive the benefits of digital technology.

These two products offer more than just a sleek design. Exceptional comfort and proven features effectively change the way people hear with an affordable hearing instrument.



User-friendly *hearing solutions*

For over five decades, Rexton has partnered with hearing health professionals to provide innovative, intelligent solutions for improved hearing.

Contact your hearing health professional to learn more about the impact the new Targa series of intermediate digital hearing instruments can have on your quality of life.



Rev 07/07 PN 0001284/b
REXTON

Targa *PLUS* & Targa *PRO*



intermediate *digital* instruments

REXTON



Our new Targa products combine advanced digital signal processing technology with superior quality. From reducing feedback to improving the sound of speech in noise, the Targa PRO and Targa PLUS reveal the knowledge, experience and confidence Rexton has demonstrated to hearing instrument users since our introduction to the industry over 50 years ago.

Directional Microphone*

With the use of Targa's directional microphone, listening in noisy environments is more comfortable than before. Both Targa PRO and Targa PLUS BTE instruments now have the capability of a directional microphone offering improvements in difficult situations where you need to discriminate speech from noise. Directional Microphones are recommended for users with a high demand for speech understanding and unpredictable listening environments.

Feedback Cancellation

With the Targa's feedback cancellation system, wearers of Targa hearing instruments are able to enjoy a more comfortable listening experience with reduced occurrence of feedback. The enhancements automatically reduce feedback and help to avoid annoying sounds.

Telephone Compatibility

Cordless and cellular phones can interfere with the performance of some hearing instruments. Targa's digital processor is equipped with a shield-protected microphone, helping to reduce acoustical interference that may occur during phone use. By changing the programmable features in your instrument, a hearing health professional can maximize the performance of Targa's telecoil option for optimum clarity while on the telephone. Consult your hearing professional about your daily telephone use. He or she will help you determine if a telecoil is appropriate for your Targa digital instrument.

Personalized Listening Situation

Targa digital hearing instruments allow hearing health professionals to program up to three personalized listening environments based upon a patient's amplification needs, personal preferences and size of the hearing instrument.

Memory Selection Alert Tone

The Memory Selection Alert Tone, standard on all Targa product offerings, indicates which personalized listening situation has been selected.

Low Battery Alert Tone

Standard on all Targa product offerings, the low battery alert tone helps to remind you when the battery needs changing. The hearing instrument will emit a series of beeps, indicating the battery is low, continuing to do so until the battery is replaced or dead.

Targa PLUS and PRO FREE

Targa PLUS and PRO FREE hearing instruments offer a comfortable, non-occluding, natural sounding hearing instrument. This micro BTE has a contemporary housing that offers maximum cosmetic appeal and provides a quick and accurate approach to instrument fitting.



Model Availability

All Targa digital circuits are available in a full range of cosmetically appealing models, from Completely-in-the-Canal (CIC) to and Open Fit Behind-the-Ear (BTE). Ask your hearing professional which one is right for you.

Standard Features

- 3 Programmable Memories
- Directional Microphone*
- Feedback Cancellation
- Telecoil*
- Power On Delay
- Low Battery Warning
- Programmable On/Off Switch

Targa PLUS Features

- 4 Compression Channels
- 4 Frequency Bands
- 4-Channel Noise Reduction

Targa PRO Features

- 2 Compression Channels
- 2 Frequency Bands
- Microphone Noise Reduction

*Standard on BTE products, optional on ITE products.



Quality digital hearing instruments for a *better quality of life.*